

Date of issue: _____

Patient:

First name:	_____
Surname:	_____
Date of birth:	_____

Practice:

Dental laboratory:

Technician: _____

Date received: _____

Date sent: _____

Date sent: _____

Implant location:

18	17	16	15	14	13	12	11		21	22	23	24	25	26	27	28
∅	∅	∅	∅	∅	∅	∅	∅	Implant dimensions ¹ (label) After implant placement the adhesive labels included in the implant packaging can be used for documentation.	∅	∅	∅	∅	∅	∅	∅	∅
L	L	L	L	L	L	L	L		L	L	L	L	L	L	L	L
								Abutment (Refer to label for description)								
								Prosthetic restoration ²								
∅	∅	∅	∅	∅	∅	∅	∅	Implant dimensions ¹ (label) After implant placement the adhesive labels included in the implant packaging can be used for documentation.	∅	∅	∅	∅	∅	∅	∅	∅
L	L	L	L	L	L	L	L		L	L	L	L	L	L	L	L
								Abutment (Refer to label for description)								
								Prosthetic restoration ²								
48	47	46	45	44	43	42	41		31	32	33	34	35	36	37	38

¹ Implant dimensions / Series of abutments (S-M-L concept):

Implant length	7.0 / 9.0 / 11.0 / 13.0 / 15.0 mm				
Implant diameter	3.3 mm	3.7 mm	4.2 mm	4.8 mm	5.5 mm
Series of abutments	S	M	L		

All components are marked with the relevant series of abutment S, M or L. The 3.3 mm and 5.5 mm ∅ implants are not available in 7.0 mm length.

² Prosthetic restoration:

cement-retained / screw-retained		removable				
Single restoration	Bridgework	Resilient bar / Bar attachment	Milled bar	Telescope	Ball attachment	tioloc
E	B	S	G	T	K	L

Connector geometry: Conical connector Platform connector

Material information:

Alloy: _____ Adhesive: _____ Ceramic: _____
 Solder: _____ Acrylic: _____ Tooth shade: _____

^{*} For tiologic® TWINFIT only.