

Date of issue: \_\_\_\_\_

**Patient:**

First name:	
Surname:	
Date of birth:	

**Practice:**

**Dental laboratory:**

Technician: \_\_\_\_\_

Date received: \_\_\_\_\_

Date sent: \_\_\_\_\_

Date sent: \_\_\_\_\_

**Implant location:**

18	17	16	15	14	13	12	11		21	22	23	24	25	26	27	28
∅	∅	∅	∅	∅	∅	∅	∅	Implant dimensions <sup>1</sup> (label) After implant placement the adhesive labels included in the implant packaging can be used for documentation.	∅	∅	∅	∅	∅	∅	∅	∅
L	L	L	L	L	L	L	L		L	L	L	L	L	L	L	L
								Abutment (Refer to label for description)								
								Prosthetic restoration <sup>2</sup>								
∅	∅	∅	∅	∅	∅	∅	∅	Implant dimensions <sup>1</sup> (label) After implant placement the adhesive labels included in the implant packaging can be used for documentation.	∅	∅	∅	∅	∅	∅	∅	∅
L	L	L	L	L	L	L	L		L	L	L	L	L	L	L	L
								Abutment (Refer to label for description)								
								Prosthetic restoration <sup>2</sup>								
48	47	46	45	44	43	42	41		31	32	33	34	35	36	37	38

**<sup>1</sup> Implant dimensions / Series of abutments (S-M-L concept):**

Implant length	7.0 / 9.0 / 11.0 / 13.0 / 15.0 mm				
Implant diameter	3.3 mm	3.7 mm	4.2 mm	4.8 mm	5.5 mm
Series of abutments	S		M		L

All components are marked with the relevant series of abutment S, M or L. The 3.3 mm and 5.5 mm ∅ implants are not available in 7.0 mm length.

**<sup>2</sup> Prosthetic restoration:**

cement-retained / screw-retained		removable				
Single restoration	Bridgework	Resilient bar / Bar attachment	Milled bar	Telescope	Ball attachment	LOCATOR®
E	B	S	G	T	K	L

Connector geometry:  Conical connector  Platform connector

**Material information:**

Alloy:  \_\_\_\_\_ Adhesive:  \_\_\_\_\_ Ceramic:  \_\_\_\_\_  
 Solder:  \_\_\_\_\_ Acrylic:  \_\_\_\_\_ Tooth shade:  \_\_\_\_\_

<sup>1</sup> For tioLogic® TWINFIT only.