

Date of issue: _____

Patient:

First name:	
Surname:	
Date of birth:	

Practice:

--

Planned prosthetic restoration:

cement-retained / screw-retained		removable				
Single restoration	Bridgework	Resilient bar / Bar attachment	Milled bar	Telescope	Ball attachment	LOCATOR®
E	B	S	G	T	K	L

Connector geometry: Conical connector Platform connector

Planned implant position / Type of implant:

18	17	16	15	14	13	12	11		21	22	23	24	25	26	27	28
∅	∅	∅	∅	∅	∅	∅	∅	Implant dimensions ¹	∅	∅	∅	∅	∅	∅	∅	∅
L	L	L	L	L	L	L	L		L	L	L	L	L	L	L	L
∅	∅	∅	∅	∅	∅	∅	∅	Implant dimensions ¹	∅	∅	∅	∅	∅	∅	∅	∅
L	L	L	L	L	L	L	L		L	L	L	L	L	L	L	L
48	47	46	45	44	43	42	41		31	32	33	34	35	36	37	38

¹ Implant dimensions / Series of abutments (S-M-L concept):

Implant length	7.0 / 9.0 / 11.0 / 13.0 / 15.0				
Implant diameter	3.3 mm	3.7 mm	4.2 mm	4.8 mm	5.5 mm
Series of abutments	S		M		L

All components are marked with the relevant series of abutment S, M or L. The 3.3 mm and 5.5 mm ∅ implants are not available in 7.0 mm length.

Pack contents practice / laboratory:

Pack contents practice:

Impressions for working models _____

Impressions for situation models _____

Photographs _____

Bite registration _____

Tooth shade _____

Pack contents dental laboratory:

Surgical stent _____

Temporary restoration _____

Quote _____

Other _____

Comments:

--

Custom tray	Stock tray	Model	Bite	Restoration	Articulator	Facebow	Photo	Finish
Upper	Upper	Upper	Upper					
Lower	Lower	Lower	Lower					

Call back request

¹ For tioLogic® TWINFIT only.