

IT'S EASY  
#CLICK&  
GO. SO  
FAST &  
RELIABLE

discovery<sup>®</sup>  
SMART sl



PRODUCT RANGE

discovery<sup>®</sup> smart sl

# GENTLE, STABLE AND EFFICIENT.

## #CLEAN & GO



discovery<sup>®</sup>  
SMART sl

## Assortments

Brackets and buccal tubes:

discovery<sup>®</sup>  
SMART sl

(15-11/21-25, 45-41/31-35)

ORTHOCAST<sup>®</sup> SL

(16/26, 46/36)

ORTHOCAST<sup>®</sup> MINI  
M-SERIES

(17/27, 47/37)

### McLaughlin-Bennett-Trevisi\* 22

Slot: 0.56 x 0.76 mm / 22 x 30



| Torque           | - 14°      | - 14°      | - 7°       | - 7°       | 0°         | + 10°      | + 17°      | + 17°      | + 10°      | 0°         | - 7°       | - 7°       | - 14°      | - 14°      |
|------------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|
| Angulation       |            |            | 0°         | 0°         | + 8°       | + 8°       | + 4°       | + 4°       | + 8°       | + 8°       | 0°         | 0°         |            |            |
| In/Out           |            |            | 1.0 mm     | 0.7 mm     | 0.7 mm     | 1.1 mm     | 0.7 mm     | 0.7 mm     | 1.1 mm     | 0.7 mm     | 0.7 mm     | 1.0 mm     |            |            |
| Width            | 2.5 mm     | 3.6 mm     | 2.8 mm     | 2.8 mm     | 2.8 mm     | 2.6 mm     | 2.8 mm     | 2.8 mm     | 2.6 mm     | 2.8 mm     | 2.8 mm     | 2.8 mm     | 3.6 mm     | 2.5 mm     |
| Distal offset    | 0°         | 10°        |            |            |            |            |            |            |            |            |            |            | 10°        | 0°         |
| REF without hook |            |            | 722-002-15 | 722-002-14 |            | 722-002-12 | 722-002-11 | 722-002-21 | 722-002-22 |            | 722-002-24 | 722-002-25 |            |            |
| REF with hook    | 822-412-17 | 822-912-16 | 722-042-15 | 722-042-14 | 722-042-13 |            |            |            |            | 722-042-23 | 722-042-24 | 722-042-25 | 822-912-26 | 822-412-27 |
| Quantity         | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  |

upper right

upper left

|   |   |   |   |   |   |   |  |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|--|---|---|---|---|---|---|---|
| 7 | 6 | 5 | 4 | 3 | 2 | 1 |  | 1 | 2 | 3 | 4 | 5 | 6 | 7 |
| 7 | 6 | 5 | 4 | 3 | 2 | 1 |  | 1 | 2 | 3 | 4 | 5 | 6 | 7 |

lower right

lower left

| Torque           | - 10°      | - 20°      | - 17°      | - 12°      | 0°         | - 6°       | - 6°       | - 6°       | - 6°       | 0°         | - 12°      | - 17°      | - 20°      | - 10°      |
|------------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|
| Angulation       |            |            | + 2°       | + 2°       | + 3°       | 0°         | 0°         | 0°         | 0°         | + 3°       | + 2°       | + 2°       |            |            |
| In/Out           |            |            | 0.7 mm     | 0.7 mm     | 0.7 mm     | 1.1 mm     | 1.1 mm     | 1.1 mm     | 1.1 mm     | 0.7 mm     | 0.7 mm     | 0.7 mm     |            |            |
| Width            | 2.5 mm     | 3.6 mm     | 2.8 mm     | 2.8 mm     | 2.8 mm     | 2.6 mm     | 2.6 mm     | 2.6 mm     | 2.6 mm     | 2.8 mm     | 2.8 mm     | 2.8 mm     | 3.6 mm     | 2.5 mm     |
| Distal offset    | 0°         | 0°         |            |            |            |            |            |            |            |            |            |            | 0°         | 0°         |
| REF without hook |            |            | 722-002-45 | 722-002-44 |            | 722-002-42 | 722-002-41 | 722-002-31 | 722-002-32 |            | 722-002-34 | 722-002-35 |            |            |
| REF with hook    | 822-412-47 | 822-912-46 | 722-042-45 | 722-042-44 | 722-042-43 |            |            |            |            | 722-042-33 | 722-042-34 | 722-042-35 | 822-912-36 | 822-412-37 |
| Quantity         | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  | 10 pieces  |

### Assortments:

|          |                 |                                                    |                                |                           |                |
|----------|-----------------|----------------------------------------------------|--------------------------------|---------------------------|----------------|
| 1 case   | upper/lower 5-5 | Hook on cuspid                                     | 20 brackets                    | in blister                | REF 722-002-01 |
| 1 case   | upper/lower 5-5 | Hook on cuspid, first bicuspid and second bicuspid | 20 brackets                    | in blister                | REF 722-042-01 |
| 1 case   | upper/lower 7-7 | Hook on cuspid                                     | 20 brackets + 8 buccal tubes   | in blister                | REF 722-002-07 |
| 1 case   | upper/lower 7-7 | Hook on cuspid, first bicuspid and second bicuspid | 20 brackets + 8 buccal tubes   | in blister                | REF 722-042-07 |
| 10 cases | upper/lower 5-5 | Hook on cuspid                                     | 200 brackets                   | in assortment storage box | REF 722-002-10 |
| 10 cases | upper/lower 5-5 | Hook on cuspid, first bicuspid and second bicuspid | 200 brackets                   | in assortment storage box | REF 722-042-10 |
| 10 cases | upper/lower 6-6 | Hook on cuspid                                     | 200 brackets + 40 buccal tubes | in assortment storage box | REF 722-002-17 |
| 10 cases | upper/lower 6-6 | Hook on cuspid, first bicuspid and second bicuspid | 200 brackets + 40 buccal tubes | in assortment storage box | REF 722-042-17 |

### Accessories:

|                                                   |                |
|---------------------------------------------------|----------------|
| sl-Instrument for discovery <sup>®</sup> smart sl | REF 023-274-00 |
| discovery <sup>®</sup> smart sl demo model        | REF 892-000-36 |

Opening instrument



D  
DENTAURUM  
1886

\* The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin-Bennett and Trevisi.