

## Test Certificate

---

Test Cytotoxicity, L 929-Proliferation  
DIN EN ISO 10993-5, -12, ISO 9363-1, LM SOP 4-06-01

Test material Cobalt chromium alloy remanium 2001, material: Co 63 %, Cr 23 %, Mo 7 %, W 4,3 %, order no. 102-600-01


Manufacturer Dentaurem J.P. Winkelstroeter KG, Turnstrasse 31, D-75228 Ispringen

Test number 011516-20 (the complete report is available at the client)

Test result The present results indicate that the cobalt chromium alloy remanium 2001 does not release cytotoxic substances.



Dr. Madlon Timme  
Study Director



Dr. Stephan Rossberger