

BIOCOMPATIBILITY CERTIFICATE

Testmaterial: REMATITAN[®] cast metal
Titan grade 1 according to DIN 17850/ISO 5832-2/ASTM F67

Supplier: DENTAURUM J.P. Winkelstroeter KG
Turnstraße 31 D-75228 Ispringen

Study performed:

The following study was performed in order to determine the biocompatibility of the dental material according to ISO 10993-5, 1992: „Test on cytotoxicity; *in vitro* methods“. The material was produced at DENTAURUM.

CYTOTOXICITY

The cytotoxic effect of the test material was examined in a cell culture test *in vitro* using L 929-cells (fibroblast cell line). Examination of materials for medical devices with the extraction method and protein determination.

Result: REMATITAN[®] did not show any cytotoxic effect.

anawa BIOSERVICE Scientific Laboratories GmbH Munich


Behringstraße 6

D-82152 Planegg



Dr. Albrecht Poth

Date: April 26, 1996



Dr. Naveed Honarvar

Date: April 26, 1996